

## ClariCELL™ BTK [T474M] Kinase Assay Service

### Description

The ClariCELL™ BTK [T474M] Kinase Assay quantifies auto-phosphorylation of human full-length BTK [T474M], a mutant form of BTK that is less sensitive to ibrutinib, in human cells. The assay is useful to determine potencies of small-molecule inhibitors against the specified kinase in the context of a cellular environment. Compound testing services are available utilizing the assay.

### Overview

Human Embryonic Kidney (HEK 293) cells transiently expressing sequence verified human full-length BTK [T474M] are exposed to test compound or control, then lysed to release cellular proteins. BTK [T474M] is captured onto an assay plate, and the extent of autophosphorylation is quantified by ELISA using an antibody specific for the phosphorylation event. Cells expressing kinase deficient BTK [K430R] are also utilized as controls to calculate the % inhibition of test compounds.

### Assay Validation

#### BTK [T474M] Expression in Cells




Figure 1: BTK [T474M] or kinase dead (K430R) BTK was expressed transiently in 293 cells. Following cell lysis, an IP Western was performed with appropriate antibodies to capture and detect the total mutant BTK protein.

#### BTK [T474M] Autophosphorylation in Cells




Figure 2: BTK [T474M] or kinase dead (K430R) BTK was expressed transiently in 293 cells. Following cell lysis, an IP Western was performed with appropriate antibodies to capture and detect phospho-BTK [T474M] protein.

#### Quantification of Phosphorylation




Figure 3: BTK [T474M] or kinase dead (K430R) BTK was expressed transiently in 293 cells. Following cell lysis, an ELISA was performed to quantify the extent of auto-phosphorylation of BTK [T474M].

#### Reference Inhibitor Data




Figure 4: An autophosphorylation assay was performed for BTK and BTK [T474M] in the presence of ibrutinib, a BTK inhibitor. % inhibition data were plotted to determine EC<sub>50</sub>s.